# Two Distinct Activation Methods Yield Clinical-Scale Expansion of Peripheral Blood Derived $\gamma\delta$ T Cells

• Xuri™:

Medium usage:



### **Authors:**

Chengkang Zhang, Amy Burkall, Elizabeth Kagan, Karen Roderick, Erin Brooks, Kenneth Low, Alexis Bossie. Lonza Inc., 9900 Medical Center Drive, Rockville, MD 20850, United States

## **Abstract**

Gamma delta ( $\gamma\delta$ ) T cells have inherent ability to infiltrate solid tumors and can directly recognize and kill transformed cells independently of HLA-antigen presentation. Moreover,  $\gamma\delta$  T cells do not cause graft-versushost disease and provide a promising platform for the development of T-cell therapies targeting solid tumors. However, due to the low prevalence of  $\gamma\delta$  T cells in peripheral blood, it remains a challenge to generate sufficient numbers of  $\gamma\delta$  T cells *ex vivo* to produce a clinical dose.

In this study, we demonstrate that billions of  $\gamma\delta$  T cells can be generated from peripheral blood mononuclear cells (PBMCs) with two distinct activation methods, with TheraPEAK® T-VIVO® Medium – a chemically defined, non-animal origin, serum-free medium, in a rocking motion platform. The first method uses zoledronic acid as the activating agent and more than 1 x 10 $^{9}$  V $\delta$ 2+ T cells are produced within 14 days starting from 0.5 x 10<sup>9</sup> PBMCs. In the second approach either  $\alpha\beta$  T cells are specifically depleted from PBMCs or  $\gamma\delta$  T cells are isolated via negative selection prior to anti-CD3 and anti-CD28 co-stimulation, and billions of  $\gamma\delta$  T cells – including both V $\delta$ 1+ and V $\delta$ 2+ T cell subsets – are produced from PBMCs source. The expanded  $\gamma\delta$  T cells exhibit innate cytotoxicity towards the K562 cell line and produce cytokines including IFN-γ and TNF- $\alpha$  in vitro. Both activation methods and expansion protocols may be utilized to generate enough  $\gamma\delta$  T cells to support clinical applications.

**Summary:** In this study, we first demonstrate robust expansion of billions of  $\gamma\delta$  T cells (V $\delta$ 2+) from cryopreserved PBMCs with zoledronic acid using the TheraPEAK® T-VIVO® Medium supplemented with low amount (100 IU/mL) of IL-2, without the use of human serum. Moreover, TheraPEAK® T-VIVO® Medium supports the expansion of  $\gamma\delta$  T cells (including V $\delta$ 1+ and  $V\delta 2+$  T-cell subsets) by anti-CD3/anti-CD28 co-stimulation, following simple isolation process that depletes  $\alpha\beta$  T cells and other cells from PBMCs. Importantly, the anti-CD3/anti-CD28 approach supports high lentivirus transduction efficiency. Collectively, the data indicate that TheraPEAK® T-VIVO<sup>®</sup> Medium benefits expansion of  $\gamma\delta$  T cells for the development of T-cell therapies.

## **Download Poster:**



All TheraPEAK® Products are produced according to applicable GMP standards and follow the USP/EP guidance for cell and gene therapy raw materials. It is the end user's responsibility to ensure full compliance with all regulations based on their use of Lonza's products in their specific process. All trademarks belong to Lonza, and are registered in the USA, EU and/or CH or belong to third-party owners and are used only for informational purposes. All third-party copyrights have been reproduced with permission from their owners. The user bears the sole responsibility for determining the existence of any third-party rights, as well as obtaining any necessary licenses and approvals. The information contained herein is believed to be correct. However, no warranty is made, either expressed or implied. For more details: www.lonza.com/legal.

©2024 Lonza. All rights reserved.

Vδ2+ T-Cell Expansion Using Zoledronic Acid




Figure 1. Workflow of expanding  $V\delta 2+T$  cells out of PBMCs using zoledronic acid. Cryopreserved PBMCs are thawed and activated with 6 µM zoledronic acid in T-225 flasks. On day 7, 500 x 10<sup>6</sup> viable cells are inoculated into Xuri™ 2 L Cellbag™ (Cytiva) for further expansion until day 16.


Add IL-2 (100 IU/mL) to TheraPEAK® T-VIVO® Medium in 1 L

bag (BP08-970Y) or 5 L bag (BP08-970H), store at room

#### Expand $\gamma\delta$ T Cells Using anti-CD3/anti-CD28 Co-Stimulation

tempature (protect from light).

1L bottle + 1L bag + 5L bag = 7L



**Figure 4.** Workflow of expanding  $\gamma\delta$  T cells using anti-CD3/anti-CD28 co-stimulation. Depleting  $\alpha\beta$  T cells from PBMCs or isolating  $\gamma\delta$  T cells out of PBMCs permits the expansion of  $\gamma\delta$  T cells, including both V $\delta$ 1+ and V $\delta$ 2+ T cells. Both CD3+CD56- and CD3+CD56+ cells are included for the assessment of pan- $\alpha\beta$ TCR, pan- $\gamma\delta$ TCR, V $\delta$ 1TCR and V $\delta$ 2TCR expression.

# Generating Billions of V $\delta$ 2+ T Cells from PBMCs Using the TheraPEAK® T-VIVO® Medium


		Cell viability	Total viable cells	% of total cells			% of CD3+ T cells					
	Day			CD3- CD56+	CD3+ CD56-	CD3+ CD56+	CD3- CD56-	αβΤ %	γδ <b>Τ %</b>	Vδ <b>2</b> +%	Vδ1+%	Total γδ T cells
	0	89.6%	497 x 10 <sup>6</sup>	6.1%	51.8%	0.7%	41.4%	97.6%	0.5%	0.12%	0.4%	1.3 x 10 <sup>6</sup>
Donor #1	7	74.3%	519 x 10 <sup>6</sup>	33.2%	42.0%	11.8%	13.1%	20.6%	77.2%	69.3%	4.9%	216 x 10 <sup>6</sup>
Donor #1	10	79.5%	2,840 x 10 <sup>6</sup>	43.0%	28.3%	26.2%	2.5%	7.16%	92.0%	81.6%	7.1%	1,424 x 10 <sup>6</sup>
_	16	76.0%	8,530 x 10 <sup>6</sup>	39.4%	22.9%	36.8%	0.9%	11.1%	87.5%	76.0%	10.1%	4,456 x 10 <sup>6</sup>
	0	94.2%	719 x 10 <sup>6</sup>	9.2%	59.3%	4.9%	26.5%	95.1%	4.7%	4.6%	0.1%	21.2 x 10 <sup>6</sup>
	7	75.6%	1,250 x 10 <sup>6</sup> (*)	5.2%	29.7%	58.7%	6.4%	2.7%	95.3%	89.4%	0.7%	1,053 x 10 <sup>6</sup>
Donor #2	10	82.1%	1,300 x 10 <sup>6</sup>	1.2%	26.2%	72.1%	0.5%	1.4%	98.2%	94.6%	0.4%	1,253 x 10 <sup>6</sup>
	14	94.4%	16,700 x 10 <sup>6</sup>	0.7%	36.3%	62.7%	0.3%	1.7%	98.1%	95.0%	0.3%	16,219 x 10 <sup>6</sup>
	16	91.9%	18,700 x 10 <sup>6</sup>	0.5%	35.3%	63.9%	0.3%	2.2%	97.6%	94.8%	0.2%	18,105 x 10 <sup>6</sup>
_	0	87.4%	535 x 10 <sup>6</sup>	14.5%	27.2%	0.8%	57.5%	87.9%	11.1%	10.8%	1.8%	16.6 x 10 <sup>6</sup>
	7	87.7%	1,340 x 10 <sup>6</sup> (*)	1.9%	78.5%	15.9%	3.7%	2.0%	97.5%	95.2%	0.2%	1,233 x 10 <sup>6</sup>
Donor #3	9	89.4%	2,540 x 10 <sup>6</sup>	1.3%	59.3%	39.2%	0.2%	1.4%	98.5%	95.9%	0.2%	2,464 x 10 <sup>6</sup>
	11	96.0%	11,300 x 10 <sup>6</sup>	0.7%	66.9%	32.3%	0.2%	0.9%	99.0%	97.3%	0.1%	11,098 x 10 <sup>6</sup>
	15	94.9%	31,300 x 10 <sup>6</sup>	0.9%	61.1%	37.9%	0.1%	0.7%	99.2%	97.5%	0.1%	30,739 x 10 <sup>6</sup>

Figure 2. Expansion of V $\delta$ 2+ T cells from PBMCs with TheraPEAK® T-VIVO® Medium using zoledronic acid. For donor #2 and #3, only 500 x 10<sup>6</sup> cells are inoculated into 2 L Cellbag™ on day 7 for further expansion. Cell samples are analyzed by flow cytometry (BD FACSCelesta™). Both CD3+CD56- and CD3+CD56+ cells are included for the assessment of pan- $\alpha\beta$ TCR, pan- $\gamma\delta$ TCR, V $\delta$ 1TCR and V $\delta$ 2TCR expression.

## Expansion of Both V $\delta$ 1+ and V $\delta$ 2+ T Cells

			PBMC Donor #2	2	РВМС С	Oonor #4
Process overview		Deplete αβ T cells ι Xuri™ for up to 16 c	using anti-TCRα/β ant lays	Isolate $\gamma\delta$ T cells by negative selection; Expand in T-flask for 10 days		
Days		Day 0	Day 10	Day 16	Day 0	Day 10
Cell viability	У	90.5%	89.5%	93.7%	92.0%	88.2%
Total viable	cells	135 x 10 <sup>6</sup>	3,450 x 10 <sup>6</sup>	38,400 x 10 <sup>6</sup>	6.9 x 10 <sup>6</sup>	1,112 x 10 <sup>6</sup>
	CD3- CD56+	10.7%	0.4%	0.4%	1.0%	0.2%
% of total cells	CD3+ CD56-	4.2%	51.2%	56.6%	63.9%	63.9%
	CD3+ CD56+	0.2%	47.6%	42.2%	31.4%	35.6%
	CD3- CD56-	86.9%	0.8%	0.9%	3.7%	0.3%
	αβΤ%	2.5%	2.3%	3.6%	0.1%	0.1%
% of CD3+ T cells	γδΤ %	96.1%	88.3%	90.6%	100.0%	98.6%
	Vδ <b>2</b> +%	85.4%	57.1%	68.1%	86.5%	81.4%
	Vδ1+%	2.9%	18.6%	12.9%	9.0%	10.7%
Total γδ T c	ells	5.7 x 10 <sup>6</sup>	3,040 x 10 <sup>6</sup>	33,800 x 10 <sup>6</sup>	6.6 x 10 <sup>6</sup>	1,091 x 10 <sup>6</sup>

**Figure 5.** Expansion of  $\gamma\delta$  T cells using the anti-CD3/anti-CD28 co-stimulation approach. For donor #2,  $\alpha\beta$  T cells are depleted from PBMCs by anti-TCR $\alpha/\beta$  antibodies. The remaining cells are then activated with T Cell TransAct<sup>TM</sup> (Miltenyi) and expanded for up to 16 days in Xuri<sup>TM</sup> 2 L Cellbag<sup>TM</sup>. For donor #4,  $\gamma\delta$  T cells are isolated by negative selection before activation with T Cell TransAct™ and expanded in T-flasks for 10 days. Cell samples are analyzed by flow cytometry (BD FACSCelesta™) at indicated days. Both CD3+CD56- and CD3+CD56+ cells are included for the assessment of pan- $\alpha\beta$ TCR, pan- $\gamma\delta$ TCR, V $\delta$ 1TCR and V $\delta$ 2TCR expression.






Figure 3. Cytotoxicity of the expanded T cells towards K562 cells are assessed after 24-hour co-cultivation at various Effector: Target ratios with the K562-luc2 cell line (ATCC). The killing of target K562 cells is measured by luciferase activity using the One-Glo™ Luciferase Assay System (Promega).

#### Lentivirus Transduction in $\gamma\delta$ T Cells

Starting Materials Activation Reagent		PBMCs	Isolated $\gamma\delta$ T cells	Isolated αβ T cells TransAct™	
		Zoledronic acid	TransAct™		
Transduction o	n day	8 1 1		1	
	CD3+%	89.9%	99.5%	99.5%	
Day 13	αβΤCR+%	9.3%	0.3%	99.3%	
Flow	γδΤСR+%	90.7%	99.3%	0.7%	
Cytometry	Vδ1 TCR+%	0.6%	55.6%	0%	
	Vδ2 TCR+%	91.7%	21.6%	0.4%	




Figure 6. T cells are transduced with a GFP lentivirus under different conditions. For PBMCs activated with zoledronic acid, the cells are transduced with the GFP lentivirus on day 8 when  $\gamma\delta$  T cells become majority based on flow cytometry, at Low (10) or High (40) MOIs. For MACS-isolated  $\gamma\delta$  T cells or  $\alpha\beta$  T cells that are activated using T Cell TransAct<sup>TM</sup>, the lentivirus is added on day 1 after activation at Low (10) or High (100) MOIs. On day 13, the cells are analyzed for lentivirus transduction efficiency based on GFP expression and the cell composition using antibodies targeting CD3, CD56, pan- $\alpha\beta$ TCR, pan- $\gamma\delta$ TCR, V $\delta$ 1TCR and V $\delta$ 2TCR.