

CRISPR Editing Human Primary Resting CD4+ T Cells with RNPs using Nucleofector® Technology

Developed by researchers at the University of Munich in collaboration with Synthego and Lonza.

Introduction

This protocol describes how to deliver ribonucleoprotein (RNP) complexes that consist of purified SpCas9 nuclease complexed with Synthego's chemically modified synthetic single guide RNA (sgRNA) to human primary resting CD4+ T cells. RNP delivery is accomplished using the Lonza 4D-Nucleofector® System.

Chemically modified sgRNAs are designed to resist degradation by exonucleases and prevent innate intracellular immune cascades that can lead to cell death. The present protocol by Albanese et al. may be used for the generation of single or multiple gene knockouts in human primary resting CD4+ T cells¹. This protocol was successfully used for the knock-in of GFP into different loci and may also be used for the knock-in of other genes of interest.

To learn more about the success of this protocol, please read the [associated case study](#).

Abbreviations:

CRISPR: clustered regularly interspaced short palindromic repeats
FBS: fetal bovine serum
GFP: green fluorescent protein
HDR: Homology-directed repair
ICE: Inference of CRISPR Edits
PBS: phosphate-buffered saline
PCR: polymerase chain reaction
RNP: ribonucleoprotein
sgRNA: single guide RNA
SpCas9: *Streptococcus pyogenes* CRISPR-associated protein 9
TE: Tris-EDTA

Materials Required

Material	Ordering Information
sgRNA Synthetic sgRNA Kit (Synthego Standard Modifications)	(Synthego), available here
PCR & sequencing primers	Multiple vendors
Nuclease-free 1X TE buffer	Multiple vendors (e.g., Thermo Fisher Scientific, Catalog #AM9849)
Nuclease-free water	Multiple vendors (e.g., Thermo Fisher Scientific)
1X PBS, cell culture grade	Multiple vendors (e.g., Thermo Fisher Scientific)
4D-Nucleofector® Core Unit and X Unit	Lonza, 4D-Nucleofector® Core Unit Lonza (AAF-1003B) Lonza, 4D-Nucleofector® X Unit Lonza (AAF-1003X)
P3 Primary Cell 4D-Nucleofector® X Kit S (for 20 µl transfection volume)	(Lonza, P3 Primary Cell 4D-Nucleofector® X Kit S Lonza (V4XP-3032) for 16x 20 µl transfection volume)
RPMI 1640 with GlutaMAX Supplement	Gibco RPMI 1640 Medium, GlutaMAX™ Supplement (thermofisher.com)
10% (v/v) FBS	(Sigma)
Penicillin/streptomycin (100 IU/ml)	(Thermo Fisher Scientific)
EasySep™ Human CD4+ T Cell Enrichment Kit	(Stem Cell; 19052)
IL-7 and IL-15 (2 ng/ml each)	(Peprotech)
SpCas9 Nuclease	Multiple vendors (IDT; 1081059)
Genomic Isolation Kit	Multiple vendors (e.g., Thermo Fisher Scientific)
CaCl ₂ (1mM)	Multiple vendors (e.g., Thermo Fisher Scientific)
MgCl ₂ (3mM)	Multiple vendors (e.g., Thermo Fisher Scientific)
EDTA (1mM)	Multiple vendors (e.g., Thermo Fisher Scientific)
Triton (X 100)	Multiple vendors (e.g., Thermo Fisher Scientific)
Tris (10 mM, pH 7.5)	Multiple vendors (e.g., Thermo Fisher Scientific)
Proteinase K	Multiple vendors (e.g., Thermo Fisher Scientific)
96-well flat-bottom plates	Multiple vendors (tissue culture grade)
48-well plates	Multiple vendors (tissue culture grade)
dsDNA HDR Template (optional)	Multiple vendors (e.g., Twist Bioscience)


Note: All protocols outlined have been validated using materials mentioned in this manual. Materials other than the ones outlined in our manual may require additional optimization by the user.

General Guidelines

- Wearing gloves and using nuclease-free tubes and reagents is recommended to avoid RNase contamination.
- Always maintain sterile technique, and use sterile filter pipette tips.
- All Synthego and Lonza reagents should be stored according to the manufacturer's recommendations.
- Synthetic sgRNA should be dissolved in TE buffer and diluted to a working concentration using nuclease-free water.
- RNP complexes are stable at room temperature for up to 1 hour.
- This protocol uses 16-well Nucleocuvette® Strips for the 4D-Nucleofector® X Unit. Please note that Nucleofection® Conditions are transferable between different formats.

For further information regarding the above-described guidelines please check the [Tips & Tricks](#) section at the end of this protocol.

Workflow Graphic



The above figure was adapted from the original article from Albanese et al, Nat Methods 19, 81–89 (2022) and is licensed under the [Creative Commons Attribution License 4.0](#) which permits the use, sharing, distribution, and reproduction in any format so long as appropriate credit is given to the original authors.

Protocol

1. Isolation of Primary T cells

- Heparinized blood retained in leukocyte reduction chambers from healthy donors was used to isolate CD4⁺ T cells.
- Blood cells were diluted with non-filter sterilized PBS; isolation of CD4⁺ T cells was done using the EasySep™ Rosette Human CD4⁺ T cell enrichment kits according to the manufacturer's protocols.
- For further information please check the [Tips & Tricks](#) section at the end of this protocol.

2. Pre Nucleofection®

- A. Please make sure that the Nucleofector® Solution supplement is added to the P3 Nucleofector® Solution (16.4 µl P3 Solution + 3.6 µl supplement per sample).
- B. Start 4D-Nucleofector® System and create experimental parameter file (for details see device manual).
- C. Select the Nucleofector® Program EH-100. For other program recommendations, e.g., to achieve higher efficiency or higher viability, please see the [Tips & Tricks](#) at the end of this protocol.
- D. Pre-warm an aliquot of culture medium RPMI 1640 GlutaMAX (without supplements) to 37°C (100 µl per sample); see Step 4 for recovery post Nucleofection®.
- E. Prepare cell culture plates by filling the appropriate number of wells with the desired volume of recommended media.

3. Nucleofection® Process

3.1. Dissolve and Dilute your sgRNA

- A. Briefly centrifuge your tubes containing the sgRNA to ensure that the dried RNA pellet is collected at the bottom.
- B. If you are working with a 1.5 nmol sgRNA, rehydrate sgRNA in 15 µl nuclease-free buffer (1X TE buffer) and pulse vortex for 30 seconds to ensure complete mixing. This will make a stock solution of 100 µM (100 pmol/µl) of sgRNA.
Note: If sgRNA is not fully dissolved, let it sit overnight and up to 72 hours at 4°C. If the sgRNA is not being used immediately, dissolved sgRNA should be aliquoted into 6 µl per tube and stored at -20°C. Under these conditions, the sgRNA is stable for up to 3 years if not repeatedly thawed. For quantification guidelines and further information, please check the [Tips & Tricks](#) section at the end of this protocol.
- C. The sgRNA 100 µM stock solution is ready to use. If you wish to dilute your sgRNA to a different working concentration, please use nuclease-free water, pulse vortex the sgRNA/nuclease-free water mix for 30 seconds and incubate at room temperature for 5 minutes to dissolve the sgRNA.

Option A: Knockout Generation in Resting Human CD4+ T Cells

This protocol describes amounts for one sample; multiple samples can be processed simultaneously if desired.

3.2. Assemble RNP Complex (2.5:1 sgRNA to SpCas9 ratio) and Transfect Cells

- A. Prepare the RNP by incubating synthetic sgRNAs together with SpCas9 at a ratio of 2.5:1 to form the RNP complex as recommended in the table below.

Note: SpCas9 protein is provided in solution and can be used directly. Please use sterile filtered (0.22 µm) PBS to dilute the RNP to a final concentration of 20 µM.

RNP Preparation: Components, Concentration & Volume

Component	Quantity (pmol)	Concentration (µM)	Volume (µl)
sgRNA	250	100	2.5
SpCas9*	100	62	1.6
PBS (filter sterilized)	N/A	N/A	0.9
Total Volume	N/A	N/A	5

*Please adjust volumes if a different SpCas9 concentration is used. The original protocol uses SpCas9 from IDT, alternatively, Synthego's 2NLS SpCas9 can also be used, available [here](#).

- B. Incubate RNP mix at room temperature for 10 minutes.
- C. Count an aliquot of the cells and determine cell density.
- D. Wash cells twice with non-filtered sterilized PBS and centrifuge 2×10^6 cells at 200 x g for 10 minutes at room temperature. Remove supernatant completely and resuspend the pellet carefully in 20 µl room temperature P3 Nucleofector® Solution.
- E. Mix the RNP solution with the cell suspension as described in the table below depending on your editing method.

Nucleofection® Sample Preparation for Different Editing Methods

Method	Volume of RNP Solution	Cell Resuspension Volume in P3 Nucleofector® Solution
Single sgRNA editing	5 µl of RNP	20 µl of cells
Multi-guide sgRNA editing	2 µl of each RNP complexes for each sgRNA	20 µl of cells
Co-editing of genes	0.5 µl of each RNP complexes for each sgRNA	20 µl of cells

- F. Transfer 20 µl of the cells mixed with RNPs into a well of the 16-well Nucleocuvette® Strip.
Note: For example, for single sgRNA editing, mix RNP 5 µl + cells 20 µl very well and only use 20 µl of sample for transfection.
- G. Transfect cells using Nucleofector® Program EH-100.
- H. Recover the cells according to Step 4.

Option B: Knock-in Generation in Resting Human CD4+ T Cells

This protocol describes amounts for one sample: multiple samples can be processed simultaneously if desired.

For the knock-in of GFP or any other gene of interest, a double-stranded DNA HDR donor template with around 550 bp homology arms in front 5' and 550bp after 3' the sgRNA cutting site is recommended.

3.3. In vitro digestion of knock-in HDR DNA templates to confirm no cleavage of the HDR template by the RNP (Optional)

- Mix the RNP complex (1 μ M), the HDR DNA template (1 μ g) and single-cutter restriction enzyme for the specific plasmid (positive control of DNA cleavage).
- Incubate samples at 37°C for 2 hours.
- Add 1 μ l of proteinase K (20 mg/ml).
- Incubate at 56°C for 10 minutes.
- Visualize sample on a 1% agarose gel.

3.4. Assemble RNP Complex (2.5:1 sgRNA to SpCas9 Ratio) and Transfect Cells

- Prepare the RNP by incubating synthetic sgRNAs together with SpCas9 at a ratio of 2.5:1 to form the RNP complex as recommended in the table below.

Note: SpCas9 protein is provided in solution and can be used directly. Please use sterile filtered (0.22 μ m) PBS to dilute the RNP to a final concentration of 20 μ M.

RNP and HDR Donor Preparation: Components, Concentration & Volume

Component	Quantity (pmol)	Concentration (μ M)	Volume (μ l)
sgRNA	250	100	2.5
SpCas9 ^{1*}	100	62	1.6
PBS (filter sterilized)	N/A	N/A	0.9
Total Volume	N/A	N/A	5 ^{1*}
DNA HDR donor template	0.8–1.12 pmol ^{2*} To be added after preparation of cell suspension		

^{1*} Please adjust volumes if a different SpCas9 concentration is used. The original protocol uses SpCas9 from IDT, alternatively, Synthego's 2NLS SpCas9 can also be used, available [here](#).

^{2*} The pmol will vary with the DNA HDR donor template used. Usual templates span from 1900 bp to 1350 bp which correspond to the pmol above.

- Count an aliquot of the cells and determine cell density.
- Wash cells twice with non-filter sterilized PBS, centrifuge 2×10^6 cells at 200 x g for 10 minutes at room temperature. Remove supernatant completely and resuspend the pellet carefully in 20 μ l room temperature P3 Nucleofector[®] Solution.
- Mix 5 μ l of the RNP template solution with 20 μ l of the cell suspension.

- E. Add donor DNA template (0.8-1.12 pmol).
- F. Mix well and transfer 20 µl of the sample into a well of the 16-well Nucleocuvette® Strip.
- G. Transfect cells using Nucleofector® Program EH-100.
- H. Recover the cells according to Step 4.

4. Post Nucleofection®

- A. Add 100 µl of pre-warmed RPMI 1640 GlutaMAX (without supplements) to each well in the 48-well plate.
- B. Transfer cells to 48-well plates and allow them to recover for 10 minutes at 37 °C.
- C. After recovery, cultivate and continue regular media changes until analysis as described in Step 5 below.

5. Cultivation of Primary Human CD4+ T Cells

- A. Cultivate cells in 96-well flat-bottom plates at a cell density of 1×10^6 cells/ ml.
- B. Add complete culture medium; RPMI 1640 GlutaMAX supplemented with 10% (v/v) FBS, penicillin/streptomycin (100 IU/ ml) and IL-7 and IL-15 (2 ng/ ml each) supplements.
- C. Replace complete medium (including FBS, antibiotics and IL-7, IL-15 supplements) every 3 days.

6. Genotyping Analysis

Editing efficiencies can be assessed via Sanger sequencing through a variety of analysis tools. Synthego's ICE analysis tool can be used to assess both CRISPR knockout and knock-in edits.

To assess editing efficiency with ICE:

- A. Design PCR primers compatible with [Inference of CRISPR Edits \(ICE\) analysis](#).
- B. 1 week post Nucleofection®, collect cells (5×10^4) and lyse in lysis buffer (20 µl) (1 mM CaCl₂, 3 mM MgCl₂, 1 mM EDTA, 1% Triton X 100 and 10 mM Tris, pH 7.5) with the addition of proteinase K (20 µg/ml; Thermo Fisher Scientific). Incubate the cell lysate at 65 °C for 20 minutes, followed by 95 °C for 15 minutes and store at -20°C until the PCR specific for the CRISPR/Cas9 target sites is performed.
- C. Use 1 µl of cell lysate as a PCR template and Sanger-sequence the PCR amplicons.
- D. Alternatively, you can follow the Genotyping guidelines provided in our [Genotyping protocol](#) for DNA isolation, PCR and Sanger sequencing.
- E. Conduct [Inference of CRISPR Edits \(ICE\)](#) analysis on Sanger sequences to determine editing efficiency. See our [ICE Knockout Analysis protocol](#) or [ICE Knock-in Analysis protocol](#) for instructions on how to use ICE and interpret results.

Note: ICE should be able to handle any deletion sizes and work with any donor template size. However, Synthego has only validated ICE analysis of deletions up to 40 bases for sgRNA and up to 150 bases for multi-guide experiments. For knock-in experiments, we have tested insertions of up to 270 bases long.

The larger the deletion, the shorter the available Sanger sequence that will be used for the alignment with the control, hence the higher chance that ICE will display warning messages which are indicated when the analysis displays “Succeeded” in orange with the warning message below it. The Sanger sequence quality/ length will dictate if the analysis is successful.

Alternatively, if ICE analysis is not successful in genotyping your large knock-in edit please refer to loci-specific Sanger sequencing or functional assays to determine large knock-in success.

Tips & Tricks

• Cell preparation/ cell source


Alternatively to isolating primary T cells, cryopreserved PBMCs can be used to isolate CD4+ cells or cryopreserved CD4+ T cells can be used:

PBMCs (different fill sizes), cryopreserved	CC-2702 (50M), CC-2703 (100M), CC-2704 (10M), CC-2705 (25M)
CD4+ T cells, cryopreserved	2W-200

• Nucleofector® Program

The above protocol describes the use of pulse code EH-100 to transfect resting T-cells with Nucleofector® Technology;

Depending on your desired outcome (e.g. editing efficiency, number of viable transfected cells) you might want to test additional pulses. Please test (some of) the following pulse conditions along with a no pulse control:



The P3 Primary Cell 4D-Nucleofector® X Kit S (Lonza) is supplied with the pmaxGFP™ Vector as a positive control for the Nucleofection® Process. We recommend using 0.4 µg pmaxGFP™ Vector for a transfection volume of 20 µl.

• Viability

The centrifugation step prior to Nucleofection® Process highly influences the viability of transfected cells. Reducing centrifugation speed to 90 x g can improve cell viability.

To achieve maximum viability, we highly recommend adding pre-warmed media to each well of the Nucleocuvette® Strip directly after Nucleofection® and incubate for 10 minutes at room temperature. After that, carefully transfer the cells to a prepared culture dish in appropriate seeding density.

Medium change 6 hours post Nucleofection® may lead to increased viability after transfection. Spin the culture dishes for 8 minutes at 140 x g in a tissue culture centrifuge. Carefully remove the medium and add fresh pre-equilibrated culture medium as described in the protocol.

• Cell Numbers/ Throughput

Nucleofection® Conditions are transferable between different formats. The 4D-Nucleofector® X Unit offers the possibility to transfect in 20 µl format (described in this protocol) or in 100 µl format. For upscaling to 100 µl please adjust cell number and amount of RNP as shown in the table below:

	20 µl Volume	100 µl Volume
Kit #	P3 Primary Cell 4D-Nucleofector® X Kit S (V4XP-3032 for 16 x 20 µl reaction)	P3 Primary Cell 4D-Nucleofector® X Kit L (V4XP-3012 for 12 x 100 µl reaction; V4XP-3024 for 24 x 100 µl reaction)
Cell number	1-2 x 10 ⁶	0.5-1 x 10 ⁷
Nucleofector® Solution	20 µl (P3 + supplement)	100 µl (P3 + supplement)
SpCas9	100 pmol	500 pmol
sgRNAs	40 pmol	200 pmol
HDR template	0.8-1.12 pmol	4 - 5.6 pmol

In addition, it is possible to perform transfections in 20 µl volume in 96-well format using the [4D-Nucleofector® 96-well Unit](#) or in 384-well format using the [384-well Nucleofector® System](#).

• Genome Editing in Other Cell Types

As a starting point, this protocol may also be used for activated T cells.

If you want to perform CRISPR Cas9-mediated genome editing or use a different editing system in other cell types (primary cells or cell lines), please check Lonza's [Knowledge Center](#) for detailed protocols and pulse recommendations or contact our [Scientific Support Team](#).

Additionally, Synthego has other optimized RNP protocols available using synthetic sgRNAs for other cell types such as the [Immortalized Cell Nucleofection Protocol](#) and the [iPS Cell Nucleofection Protocol](#).

• sgRNA and RNP

A. sgRNA Design

For knockouts, it is generally recommended to try 3-5 sgRNAs with unique target sequences (evaluated independently) per target that cut near (<10 bp) the edit site and check their indel percentage by [ICE](#) first. When using two sgRNAs together to create a knockout, the distance between the cutting sites must not be multiple of 3 nt. Otherwise, you remove only few amino acids and the protein can still be functional.

For knock-ins, the higher the indel efficiency of a given sgRNA, the higher the knock-in efficiency. Therefore, it is generally recommended to try 3-5 sgRNAs with unique target sequences (evaluated independently) per target that cut near (<10 bp) the edit site and check their indel percentage by [ICE](#) first.

B. sgRNA Rehydration, Storage, and Quantification

Synthetic sgRNA can also be dissolved in nuclease-free water if desired, however, we recommend using TE buffer to help stabilize pH when complexing with the nuclease. Please refer to our [sgRNA Quick Start Guide](#) and [Gene Knockout Kit v2 Quick Start Guide](#) for Rehydration, Storage, and Quantification guidelines.

C. RNP

It is important to work as quickly as possible as RNP complexes begin to become unstable after 1 hour at room temperature. Additionally, leaving cells in Nucleofector® Solution for extended periods of time may lead to reduced transfection efficiency and viability.

RNP complexes may also be stored at 4°C for up to one week, or at -20°C for up to one month. Note that RNPs stored at 4°C may become susceptible to contamination from microbial growth after long periods of time.

Dissolving RNPs in P3 Nucleofector® Solution is an option to avoid dilution of the Nucleofector® Solution. Nucleofector® Solution is animal component-free and non-toxic.

D. HDR Template

Our [Tips and Tricks guide: Design and Optimization of CRISPR Knock-in Experiments](#) offers general guidelines for planning your knock-in experiments.

The quality and the concentration of DNA HDR donor template used for Nucleofection® plays a central role for the efficiency of gene transfer. We strongly recommend the use of high-quality plasmids from vendors or, if isolated within the lab, use plasmid purification kits like EndoFree® Plasmid Kits (Qiagen, Cat. No. 12391 Giga Kit, 12362 Maxi Kit, 12381 Mega Kit). The purified DNA HDR donor template should be resuspended in deionized water or TE buffer with a concentration between 0.2-1 µg/µl. Please check the purity of each plasmid preparation by measuring the A260:A280 ratio, according to manufacturer's protocol. Other DNA donor templates could be used but will require optimization by the user.

If you wish to try alternative knock-in strategies, please be aware that we have no specific recommendations. There are limitations on the amount of donor DNA template you can use and the editing efficiency of each knock-in may vary.

Contact Technical Support

If you require further technical support, please reach out to our Scientific Support Teams:

- Synthego (support@synthego.com):
 - sgRNA and HDR template design, sgRNA resuspension
 - RNP formation
 - ICE analysis
- Lonza (Scientific.Support@lonza.com (US), Scientific.Support.EU@lonza.com (EU and International)):
 - T-cell preparation and cell handling
 - Nucleofection® Process, protocols and format

Citations

1. Albanese M et al. Rapid, efficient and activation-neutral gene editing of polyclonal primary human resting CD4+ T cells allows complex functional analyses. Nat Methods. 2022 Jan;19(1):81-89. doi: 10.1038/s41592-021-01328-8. Epub 2021 Dec 23. PMID: 34949807; PMCID: PMC8748193.

Additional Information

For an up-to-date list of all Synthego Protocols and other resources, please visit [Synthego.com/resources](https://synthego.com/resources)

For technical assistance, contact our Scientific Support Team:

Ph: 888.611.6883 Email: support@synthego.com

About Synthego

Synthego is the leading genome engineering innovation company. The company's automated, full stack genome engineering platform enables broader access to CRISPR to accelerate basic scientific discovery, uncover cures for diseases, and develop novel synthetic biology applications. Headquartered in Silicon Valley, Synthego is used by scientists from the largest global biotechnology companies and global biology universities to unlock the potential of gene editing.

Lonza Cologne GmbH – 50829 Cologne Germany

For research use only. Not for use in diagnostic procedures.

The Nucleofector® Technology is covered by patent and/or patent pending rights owned by the Lonza Group Ltd or its affiliates. All trademarks belong to Lonza, registered in USA, EU or CH or to third party owners and used only for informational purposes. All third party copyrights have been reproduced with permission or belong to Lonza. The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. The buyer assumes all risks of use and/or handling. Any user must make his own determination and satisfy himself that the products supplied by Lonza Group Ltd or its affiliates and the information and recommendations given by Lonza Group Ltd or its affiliates are (i) suitable for intended process or purpose, (ii) in compliance with environmental, health and safety regulations, and (iii) will not infringe any third party's intellectual property rights. The user bears the sole responsibility for determining the existence of any such third party rights, as well as obtaining any necessary licenses.

For more details: www.lonza.com/legal.

© 2023 Lonza. CD-OP013 03/23.